

Antibodies Product Datasheet

Product Name:	Fe65L Antibody
Catalog Number:	C14453
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APBB2. AA range:471-520
Concentration:	1 mg/ml
Gene Name:	APBB2
Protein Name:	Amyloid beta A4 precursor protein-binding family B member 2
Other Name:	APBB2; FE65L; FE65L1; Amyloid beta A4 precursor protein-binding family B member 2; Protein Fe65-like 1
Human Swiss Prot No.:	Q92870
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Fe65L Polyclonal Antibody detects endogenous levels of Fe65L protein.
Storage and Stability:	-20°C/1 year