

Antibodies Product Datasheet

Product Name:	APP-BP2 Antibody
Catalog Number:	C14456
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from APP-BP2 . at AA range: 390-470
Concentration:	1 mg/ml
Gene Name:	APPBP2
Protein Name:	Amyloid protein-binding protein 2
Other Name:	APPBP2; KIAA0228; PAT1; Amyloid protein-binding protein 2; Amyloid beta precursor protein-binding protein 2; APP-BP2; Protein interacting with APP tail 1
Human Swiss Prot No.:	Q92624
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	APP-BP2 Polyclonal Antibody detects endogenous levels of APP-BP2 protein.
Storage and Stability:	-20°C/1 year