

Antibodies Product Datasheet

Product Name:	ALS2CR4 Antibody
Catalog Number:	C14462
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ALS2CR4. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	TMEM237
Protein Name:	Transmembrane protein 237
Other Name:	TMEM237; ALS2CR4; Transmembrane protein 237; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 4 protein
Human Swiss Prot No.:	Q96Q45
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ALS2CR4 Polyclonal Antibody detects endogenous levels of ALS2CR4 protein.
Storage and Stability:	-20°C/1 year