

Antibodies Product Datasheet

Product Name:	AIG-1 Antibody
Catalog Number:	C14468
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AIG1. AA range:196-245
Concentration:	1 mg/ml
Gene Name:	AIG1
Protein Name:	Androgen-induced gene 1 protein
Other Name:	AIG1; CGI-103; Androgen-induced gene 1 protein; AIG-1
Human Swiss Prot No.:	Q9NVV5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AIG-1 Polyclonal Antibody detects endogenous levels of AIG-1 protein.
Storage and Stability:	-20°C/1 year