

Antibodies Product Datasheet

Product Name:	KANK2 Antibody
Catalog Number:	C14496
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KANK2. AA range:351-400
Concentration:	1 mg/ml
Gene Name:	KANK2
Protein Name:	KN motif and ankyrin repeat domain-containing protein 2
Other Name:	KANK2; ANKRD25; KIAA1518; MXRA3; SIP; KN motif and ankyrin repeat domain-containing protein 2; Ankyrin repeat domain-containing protein 25; Matrix-remodeling-associated protein 3; SRC-1-interacting protein; SIP; SRC-interacting protein; SRC
Human Swiss Prot No.:	Q63ZY3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	KANK2 Polyclonal Antibody detects endogenous levels of KANK2 protein.
Storage and Stability:	-20°C/1 year