

Antibodies Product Datasheet

Product Name:	ANKRD26 Antibody
Catalog Number:	C14497
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ANKRD26. AA range:791-840
Concentration:	1 mg/ml
Gene Name:	ANKRD26
Protein Name:	Ankyrin repeat domain-containing protein 26
Other Name:	ANKRD26; KIAA1074; Ankyrin repeat domain-containing protein 26
Human Swiss Prot No.:	Q9UPS8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ANKRD26 Polyclonal Antibody detects endogenous levels of ANKRD26 protein.
Storage and Stability:	-20°C/1 year