

Antibodies Product Datasheet

Product Name:	Annexin X Antibody
Catalog Number:	C14499
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ANXA10. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	ANXA10
Protein Name:	Annexin A10
Other Name:	ANXA10; ANX14; Annexin A10; Annexin-10; Annexin-14
Human Swiss Prot No.:	Q9UJ72
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Annexin X Polyclonal Antibody detects endogenous levels of Annexin X protein.
Storage and Stability:	-20°C/1 year