

Antibodies Product Datasheet

Product Name:	ApoL1 Antibody
Catalog Number:	C14539
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APOL1. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	APOL1
Protein Name:	Apolipoprotein L1
Other Name:	APOL1; APOL; Apolipoprotein L1; Apolipoprotein L; Apo-L; ApoL; Apolipoprotein L-I; ApoL-I
Human Swiss Prot No.:	O14791
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ApoL1 Polyclonal Antibody detects endogenous levels of ApoL1 protein.
Storage and Stability:	-20°C/1 year