

## Antibodies Product Datasheet

|                               |                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | ApoL3 Antibody                                                                                                        |
| <b>Catalog Number:</b>        | C14541                                                                                                                |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                     |
| <b>Applications:</b>          | IHC-p, IF(paraffin section), ELISA                                                                                    |
| <b>Dilutions:</b>             | Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.                            |
| <b>Modification:</b>          | Unmodified/Total                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from ApoL3 . at AA range: 100-180                                                         |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                               |
| <b>Gene Name:</b>             | APOL3                                                                                                                 |
| <b>Protein Name:</b>          | Apolipoprotein L3                                                                                                     |
| <b>Other Name:</b>            | APOL3; Apolipoprotein L3; Apolipoprotein L-III; ApoL-III; TNF-inducible protein CG12-1; CG12_1                        |
| <b>Human Swiss Prot No.:</b>  | O95236                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Specificity:</b>           | ApoL3 Polyclonal Antibody detects endogenous levels of ApoL3 protein.                                                 |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                          |