

Antibodies Product Datasheet

Product Name:	ApoL5 Antibody
Catalog Number:	C14543
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APOL5. AA range:384-433
Concentration:	1 mg/ml
Gene Name:	APOL5
Protein Name:	Apolipoprotein L5
Other Name:	APOL5; Apolipoprotein L5; Apolipoprotein L-V; ApoL-V
Human Swiss Prot No.:	Q9BWW9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ApoL5 Polyclonal Antibody detects endogenous levels of ApoL5 protein.
Storage and Stability:	-20°C/1 year