

Antibodies Product Datasheet

Product Name:	AQP11 Antibody
Catalog Number:	C14548
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from AQP11 . at AA range: 130-210
Concentration:	1 mg/ml
Gene Name:	AQP11
Protein Name:	Aquaporin-11
Other Name:	AQP11; AQPX1; PSEC0027; Aquaporin-11; AQP-11
Human Swiss Prot No.:	Q8NBQ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AQP11 Polyclonal Antibody detects endogenous levels of AQP11 protein.
Storage and Stability:	-20°C/1 year