

Antibodies Product Datasheet

Product Name:	AQP12 Antibody
Catalog Number:	C14549
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AQP12. AA range:231-280
Concentration:	1 mg/ml
Gene Name:	AQP12A/AQP12B
Protein Name:	Aquaporin-12A/B
Other Name:	AQP12A; AQP12; AQPX2; Aquaporin-12A; AQP-12; AQP12B; Aquaporin-12B; AQP-12B
Human Swiss Prot No.:	Q8IXF9/A6NM10
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AQP12 Polyclonal Antibody detects endogenous levels of AQP12 protein.
Storage and Stability:	-20°C/1 year