

Antibodies Product Datasheet

Product Name:	Arginase II Antibody
Catalog Number:	C14560
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARG2. AA range:305-354
Concentration:	1 mg/ml
Gene Name:	ARG2
Protein Name:	Arginase-2 mitochondrial
Other Name:	ARG2; Arginase-2; mitochondrial; Kidney-type arginase; Non-hepatic arginase; Type II arginase
Human Swiss Prot No.:	P78540
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Arginase II Polyclonal Antibody detects endogenous levels of Arginase II protein.
Storage and Stability:	-20°C/1 year