

Antibodies Product Datasheet

Product Name:	Arylsulfatase F Antibody
Catalog Number:	C14568
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARSF. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	ARSF
Protein Name:	Arylsulfatase F
Other Name:	ARSF; Arylsulfatase F; ASF
Human Swiss Prot No.:	P54793
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Arylsulfatase F Polyclonal Antibody detects endogenous levels of Arylsulfatase F protein.
Storage and Stability:	-20°C/1 year