

Antibodies Product Datasheet

Product Name:	Ataxin-2 Antibody
Catalog Number:	C14578
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATXN2. AA range:731-780
Concentration:	1 mg/ml
Gene Name:	ATXN2
Protein Name:	Ataxin-2
Other Name:	ATXN2; ATX2; SCA2; TNRC13; Ataxin-2; Spinocerebellar ataxia type 2 protein; Trinucleotide repeat-containing gene 13 protein
Human Swiss Prot No.:	Q99700
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ataxin-2 Polyclonal Antibody detects endogenous levels of Ataxin-2 protein.
Storage and Stability:	-20°C/1 year