

Antibodies Product Datasheet

Product Name:	Ataxin-2L Antibody
Catalog Number:	C14580
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATXN2L. AA range:571-620
Concentration:	1 mg/ml
Gene Name:	ATXN2L
Protein Name:	Ataxin-2-like protein
Other Name:	ATXN2L; A2D; A2LG; A2LP; A2RP; Ataxin-2-like protein; Ataxin-2 domain protein; Ataxin-2-related protein
Human Swiss Prot No.:	Q8WWM7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ataxin-2L Polyclonal Antibody detects endogenous levels of Ataxin-2L protein.
Storage and Stability:	-20°C/1 year