

Antibodies Product Datasheet

Product Name:	Ataxin-7L1 Antibody
Catalog Number:	C14583
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATXN7L1. AA range:441-490
Concentration:	1 mg/ml
Gene Name:	ATXN7L1
Protein Name:	Ataxin-7-like protein 1
Other Name:	ATXN7L1; ATXN7L4; KIAA1218; Ataxin-7-like protein 1; Ataxin-7-like protein 4
Human Swiss Prot No.:	Q9ULK2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ataxin-7L1 Polyclonal Antibody detects endogenous levels of Ataxin-7L1 protein.
Storage and Stability:	-20°C/1 year