

Antibodies Product Datasheet

Product Name:	ATP5F1 Antibody
Catalog Number:	C14590
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ATP5F1 . at AA range: 130-210
Concentration:	1 mg/ml
Gene Name:	ATP5F1
Protein Name:	ATP synthase subunit b mitochondrial
Other Name:	ATP5F1; ATP synthase subunit b; mitochondrial; ATPase subunit b
Human Swiss Prot No.:	P24539
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5F1 Polyclonal Antibody detects endogenous levels of ATP5F1 protein.
Storage and Stability:	-20°C/1 year