

Antibodies Product Datasheet

Product Name:	ATP5C1 Antibody
Catalog Number:	C14596
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATP5C1. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	ATP5C1
Protein Name:	ATP synthase subunit gamma mitochondrial
Other Name:	ATP5C1; ATP5C; ATP5CL1; ATP synthase subunit gamma; mitochondrial; F-ATPase gamma subunit
Human Swiss Prot No.:	P36542
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5C1 Polyclonal Antibody detects endogenous levels of ATP5C1 protein.
Storage and Stability:	-20°C/1 year