

Antibodies Product Datasheet

Product Name:	ATP5G2 Antibody
Catalog Number:	C14598
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATP5G2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	ATP5G2
Protein Name:	ATP synthase lipid-binding protein mitochondrial
Other Name:	ATP5G2; PSEC0033; ATP synthase lipid-binding protein; mitochondrial; ATP synthase proteolipid P2; ATPase protein 9; ATPase subunit c
Human Swiss Prot No.:	Q06055
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5G2 Polyclonal Antibody detects endogenous levels of ATP5G2 protein.
Storage and Stability:	-20°C/1 year