

Antibodies Product Datasheet

Product Name:	ATP5G3 Antibody
Catalog Number:	C14599
Reactivity:	Human, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATP5G3. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	ATP5G3
Protein Name:	ATP synthase lipid-binding protein mitochondrial
Other Name:	ATP5G3; ATP synthase lipid-binding protein; mitochondrial; ATP synthase proteolipid P3; ATPase protein 9; ATPase subunit c
Human Swiss Prot No.:	P48201
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5G3 Polyclonal Antibody detects endogenous levels of ATP5G3 protein.
Storage and Stability:	-20°C/1 year