

Antibodies Product Datasheet

Product Name:	ATP5L2 Antibody
Catalog Number:	C14604
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATP5L2. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	ATP5L2
Protein Name:	ATP synthase subunit g 2 mitochondrial
Other Name:	ATP5L2; ATP5K2; ATP synthase subunit g 2; mitochondrial; ATPase subunit g 2
Human Swiss Prot No.:	Q7Z4Y8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5L2 Polyclonal Antibody detects endogenous levels of ATP5L2 protein.
Storage and Stability:	-20°C/1 year