

Antibodies Product Datasheet

Product Name:	ATP5S Antibody
Catalog Number:	C14607
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATP5S. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	ATP5S
Protein Name:	ATP synthase subunit s mitochondrial
Other Name:	ATP5S; ATPW; ATP synthase subunit s; mitochondrial; ATP synthase-coupling factor B; FB; Mitochondrial ATP synthase regulatory component factor B
Human Swiss Prot No.:	Q99766
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP5S Polyclonal Antibody detects endogenous levels of ATP5S protein.
Storage and Stability:	-20°C/1 year