

Antibodies Product Datasheet

Product Name:	ABCB10 Antibody
Catalog Number:	C14619
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABCB10. AA range:441-490
Concentration:	1 mg/ml
Gene Name:	ABCB10
Protein Name:	ATP-binding cassette sub-family B member 10 mitochondrial
Other Name:	ABCB10; ATP-binding cassette sub-family B member 10; mitochondrial; ATP-binding cassette transporter 10; ABC transporter 10 protein; Mitochondrial ATP-binding cassette 2; M-ABC2
Human Swiss Prot No.:	Q9NRK6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABCB10 Polyclonal Antibody detects endogenous levels of ABCB10 protein.
Storage and Stability:	-20°C/1 year