

Antibodies Product Datasheet

Product Name:	ABCD4 Antibody
Catalog Number:	C14627
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABCD4. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	ABCD4
Protein Name:	ATP-binding cassette sub-family D member 4
Other Name:	ABCD4; PXMP1L; ATP-binding cassette sub-family D member 4; PMP70-related protein; P70R; Peroxisomal membrane protein 1-like; PXMP1-L; Peroxisomal membrane protein 69; PMP69
Human Swiss Prot No.:	O14678
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ABCD4 Polyclonal Antibody detects endogenous levels of ABCD4 protein.
Storage and Stability:	-20°C/1 year