

Antibodies Product Datasheet

Product Name:	DDX19B Antibody
Catalog Number:	C14648
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DDX19B. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	DDX19B
Protein Name:	ATP-dependent RNA helicase DDX19B
Other Name:	DDX19B; DBP5; DDX19; TDBP; ATP-dependent RNA helicase DDX19B; DEAD box RNA helicase DEAD5; DEAD box protein 19B
Human Swiss Prot No.:	Q9UMR2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DDX19B Polyclonal Antibody detects endogenous levels of DDX19B protein.
Storage and Stability:	-20°C/1 year