

Antibodies Product Datasheet

Product Name:	DDX54 Antibody
Catalog Number:	C14657
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from DDX54 . at AA range: 570-650
Concentration:	1 mg/ml
Gene Name:	DDX54
Protein Name:	ATP-dependent RNA helicase DDX54
Other Name:	DDX54; ATP-dependent RNA helicase DDX54; ATP-dependent RNA helicase DP97; DEAD box RNA helicase 97 kDa; DEAD box protein 54 Q8TDD1
Human Swiss Prot No.:	Q8TDD1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DDX54 Polyclonal Antibody detects endogenous levels of DDX54 protein.
Storage and Stability:	-20°C/1 year