

Antibodies Product Datasheet

Product Name:	DDX8 Antibody
Catalog Number:	C14659
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DHX8. AA range:391-440
Concentration:	1 mg/ml
Gene Name:	DHX8
Protein Name:	ATP-dependent RNA helicase DHX8
Other Name:	DHX8; DDX8; ATP-dependent RNA helicase DHX8; DEAH box protein 8; RNA helicase HRH1
Human Swiss Prot No.:	Q14562
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DDX8 Polyclonal Antibody detects endogenous levels of DDX8 protein.
Storage and Stability:	-20°C/1 year