

Antibodies Product Datasheet

Product Name:	Atg16 Antibody
Catalog Number:	C14670
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Atg16 . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	ATG16L1
Protein Name:	Autophagy-related protein 16-1
Other Name:	ATG16L1; APG16L; Autophagy-related protein 16-1; APG16-like 1
Human Swiss Prot No.:	Q676U5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Atg16 Polyclonal Antibody detects endogenous levels of Atg16 protein.
Storage and Stability:	-20°C/1 year