

Antibodies Product Datasheet

Product Name:	BAGE2 Antibody
Catalog Number:	C14675
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BAGE2. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	BAGE2
Protein Name:	B melanoma antigen 2
Other Name:	BAGE2; B melanoma antigen 2; Cancer/testis antigen 2.2; CT2.2
Human Swiss Prot No.:	Q86Y30
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BAGE2 Polyclonal Antibody detects endogenous levels of BAGE2 protein.
Storage and Stability:	-20°C/1 year