

Antibodies Product Datasheet

Product Name:	BNIP-2 Antibody
Catalog Number:	C14696
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, FCM, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BNIP2. AA range:265-314
Concentration:	1 mg/ml
Gene Name:	BNIP2
Protein Name:	BCL2/adenovirus E1B 19 kDa protein-interacting protein 2
Other Name:	BNIP2; NIP2; BCL2/adenovirus E1B 19 kDa protein-interacting protein 2
Human Swiss Prot No.:	Q12982
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BNIP-2 Polyclonal Antibody detects endogenous levels of BNIP-2 protein.
Storage and Stability:	-20°C/1 year