

Antibodies Product Datasheet

Product Name:	BCL2L12 Antibody
Catalog Number:	C14705
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from BCL2L12 . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	BCL2L12
Protein Name:	Bcl-2-like protein 12
Other Name:	BCL2L12; BPR; Bcl-2-like protein 12; Bcl2-L-12; Bcl-2-related proline-rich protein
Human Swiss Prot No.:	Q9HB09
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BCL2L12 Polyclonal Antibody detects endogenous levels of BCL2L12 protein.
Storage and Stability:	-20°C/1 year