

Antibodies Product Datasheet

Product Name:	GCNT3 Antibody
Catalog Number:	C14708
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GCNT3. AA range:226-275
Concentration:	1 mg/ml
Gene Name:	GCNT3
Protein Name:	Beta-1,3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 3
Other Name:	GCNT3; Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 3; C2GnT-mucin type; C2GnT-M; hC2GnT-M; Core 2/core 4 beta-1,6-N-acetylglucosaminyltransferase; C2/4GnT
Human Swiss Prot No.:	O95395
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GCNT3 Polyclonal Antibody detects endogenous levels of GCNT3 protein.
Storage and Stability:	-20°C/1 year