

Antibodies Product Datasheet

Product Name:	Manic Fringe Antibody
Catalog Number:	C14719
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MFNG. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	MFNG
Protein Name:	Beta-1,3-N-acetylglucosaminyltransferase manic fringe
Other Name:	MFNG; Beta-1; 3-N-acetylglucosaminyltransferase manic fringe; O-fucosylpeptide 3-beta-N-acetylglucosaminyltransferase
Human Swiss Prot No.:	O00587
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Manic Fringe Polyclonal Antibody detects endogenous levels of Manic Fringe protein.
Storage and Stability:	-20°C/1 year