

Antibodies Product Datasheet

Product Name:	Galactosidase β Antibody
Catalog Number:	C14732
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human Galactosidase β .
Concentration:	1 mg/ml
Gene Name:	GLB1
Protein Name:	Beta-galactosidase
Other Name:	GLB1; ELNR1; Beta-galactosidase; Acid beta-galactosidase; Lactase; Elastin receptor 1
Human Swiss Prot No.:	P16278
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Galactosidase β Polyclonal Antibody detects endogenous levels of Galactosidase β protein.
Storage and Stability:	-20°C/1 year