

Antibodies Product Datasheet

Product Name:	CEL Antibody
Catalog Number:	C14748
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CEL.
Concentration:	1 mg/ml
Gene Name:	CEL
Protein Name:	Bile salt-activated lipase
Other Name:	CEL; BAL; Bile salt-activated lipase; BAL; Bile salt-stimulated lipase; BSSL; Bucelipase; Carboxyl ester lipase; Cholesterol esterase; Pancreatic lysophospholipase; Sterol esterase
Human Swiss Prot No.:	P19835
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CEL Polyclonal Antibody detects endogenous levels of CEL protein.
Storage and Stability:	-20°C/1 year