

Antibodies Product Datasheet

Product Name:	Brp44 Antibody
Catalog Number:	C14780
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BRP44. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	MPC2
Protein Name:	Mitochondrial pyruvate carrier 2
Other Name:	MPC2; BRP44; Mitochondrial pyruvate carrier 2; Brain protein 44
Human Swiss Prot No.:	O95563
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Brp44 Polyclonal Antibody detects endogenous levels of Brp44 protein.
Storage and Stability:	-20°C/1 year