

Antibodies Product Datasheet

Product Name:	BAIAP2L2 Antibody
Catalog Number:	C14788
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BAIAP2L2. AA range:111-160
Concentration:	1 mg/ml
Gene Name:	BAIAP2L2
Protein Name:	Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 2
Other Name:	BAIAP2L2; Brain-specific angiogenesis inhibitor 1-associated protein 2-like protein 2; BAI1-associated protein 2-like protein 2; Planar intestinal- and kidney-specific BAR domain protein; Pinkbar
Human Swiss Prot No.:	Q6UXY1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BAIAP2L2 Polyclonal Antibody detects endogenous levels of BAIAP2L2 protein.
Storage and Stability:	-20°C/1 year