

Antibodies Product Datasheet

Product Name:	Tip5 Antibody
Catalog Number:	C14806
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BAZ2A. AA range:1281-1330
Concentration:	1 mg/ml
Gene Name:	BAZ2A
Protein Name:	Bromodomain adjacent to zinc finger domain protein 2A
Other Name:	BAZ2A; KIAA0314; TIP5; Bromodomain adjacent to zinc finger domain protein 2A; Transcription termination factor I-interacting protein 5; TTF-I-interacting protein 5; Tip5; hWALp3
Human Swiss Prot No.:	Q9UIF9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Tip5 Polyclonal Antibody detects endogenous levels of Tip5 protein.
Storage and Stability:	-20°C/1 year