

Antibodies Product Datasheet

Product Name:	CHP2 Antibody
Catalog Number:	C14811
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CHP2. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	CHP2
Protein Name:	Calcineurin B homologous protein 2
Other Name:	CHP2; HCA520; Calcineurin B homologous protein 2; Hepatocellular carcinoma-associated antigen 520
Human Swiss Prot No.:	O43745
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHP2 Polyclonal Antibody detects endogenous levels of CHP2 protein.
Storage and Stability:	-20°C/1 year