

Antibodies Product Datasheet

Product Name:	DSCR 1 Antibody
Catalog Number:	C14815
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RCAN1. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	RCAN1
Protein Name:	Calcipressin-1
Other Name:	RCAN1; ADAPT78; CSP1; DSC1; DSCR1; Calcipressin-1; Adapt78; Down syndrome critical region protein 1; Myocyte-enriched calcineurin-interacting protein 1; MCIP1; Regulator of calcineurin 1
Human Swiss Prot No.:	P53805
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DSCR 1 Polyclonal Antibody detects endogenous levels of DSCR 1 protein.
Storage and Stability:	-20°C/1 year