

Antibodies Product Datasheet

Product Name:	CaMKII δ Antibody
Catalog Number:	C14834
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-300 ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CaMKII δ .
Concentration:	1 mg/ml
Gene Name:	CAMK2D
Protein Name:	Calcium/calmodulin-dependent protein kinase type II subunit delta
Other Name:	CAMK2D; CAMKD; Calcium/calmodulin-dependent protein kinase type II subunit delta; CaM kinase II subunit delta; CaMK-II subunit delta
Human Swiss Prot No.:	Q13557
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CaMKII δ Polyclonal Antibody detects endogenous levels of CaMKII δ protein.
Storage and Stability:	-20°C/1 year