

Antibodies Product Datasheet

Product Name:	CHP Antibody
Catalog Number:	C14851
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human CHP.
Concentration:	1 mg/ml
Gene Name:	CHP1
Protein Name:	Calcineurin B homologous protein 1
Other Name:	CHP1; CHP; Calcineurin B homologous protein 1; Calcineurin B-like protein; Calcium-binding protein CHP; Calcium-binding protein p22; EF-hand calcium-binding domain-containing protein p22
Human Swiss Prot No.:	Q99653
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CHP Polyclonal Antibody detects endogenous levels of CHP protein.
Storage and Stability:	-20°C/1 year