

Antibodies Product Datasheet

Product Name:	iPLA2 γ Antibody
Catalog Number:	C14856
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PNPLA8. AA range:691-740
Concentration:	1 mg/ml
Gene Name:	PNPLA8
Protein Name:	Calcium-independent phospholipase A2-gamma
Other Name:	PNPLA8; IPLA22; IPLA2G; BM-043; Calcium-independent phospholipase A2-gamma; Intracellular membrane-associated calcium-independent phospholipase A2 gamma; iPLA2-gamma; PNPLA-gamma; Patatin-like phospholipase domain-containing protein 8; iPLA
Human Swiss Prot No.:	Q9NP80
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	iPLA2 γ Polyclonal Antibody detects endogenous levels of iPLA2 γ protein.
Storage and Stability:	-20°C/1 year