

Antibodies Product Datasheet

Product Name:	Calpain 9 Antibody
Catalog Number:	C14878
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CAPN9. AA range:481-530
Concentration:	1 mg/ml
Gene Name:	CAPN9
Protein Name:	Calpain-9
Other Name:	CAPN9; NCL4; Calpain-9; Digestive tract-specific calpain; New calpain 4; nCL-4; Protein CG36
Human Swiss Prot No.:	O14815
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Calpain 9 Polyclonal Antibody detects endogenous levels of Calpain 9 protein.
Storage and Stability:	-20°C/1 year