

Antibodies Product Datasheet

Product Name:	ARPP-21 Antibody
Catalog Number:	C14899
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ARPP21. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	ARPP21
Protein Name:	cAMP-regulated phosphoprotein 21
Other Name:	ARPP21; TARPP; cAMP-regulated phosphoprotein 21; ARPP-21; Thymocyte cAMP-regulated phosphoprotein
Human Swiss Prot No.:	Q9UBL0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARPP-21 Polyclonal Antibody detects endogenous levels of ARPP-21 protein.
Storage and Stability:	-20°C/1 year