

Antibodies Product Datasheet

Product Name:	MRP3 Antibody
Catalog Number:	C14907
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABCC3. AA range:971-1020
Concentration:	1 mg/ml
Gene Name:	ABCC3
Protein Name:	Canalicular multispecific organic anion transporter 2
Other Name:	ABCC3; CMOAT2; MLP2; MRP3; Canalicular multispecific organic anion transporter 2; ATP-binding cassette sub-family C member 3; Multi-specific organic anion transporter D; MOAT-D; Multidrug resistance-associated protein 3
Human Swiss Prot No.:	O15438
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRP3 Polyclonal Antibody detects endogenous levels of MRP3 protein.
Storage and Stability:	-20°C/1 year