

Antibodies Product Datasheet

Product Name:	CAGE-1 Antibody
Catalog Number:	C14912
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CAGE1. AA range:711-760
Concentration:	1 mg/ml
Gene Name:	CAGE1
Protein Name:	Cancer-associated gene 1 protein
Other Name:	CAGE1; CTAG3; Cancer-associated gene 1 protein; CAGE-1; Cancer/testis antigen 3; CT3
Human Swiss Prot No.:	Q8TC20
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CAGE-1 Polyclonal Antibody detects endogenous levels of CAGE-1 protein.
Storage and Stability:	-20°C/1 year