

Antibodies Product Datasheet

Product Name:	CLIP-170 Antibody
Catalog Number:	C14915
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLIP1. AA range:1291-1340
Concentration:	1 mg/ml
Gene Name:	CLIP1
Protein Name:	CAP-Gly domain-containing linker protein 1
Other Name:	CLIP1; CYLN1; RSN; CAP-Gly domain-containing linker protein 1; Cytoplasmic linker protein 1; Cytoplasmic linker protein 170 alpha-2; CLIP-170; Reed-Sternberg intermediate filament-associated protein; Restin
Human Swiss Prot No.:	P30622
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLIP-170 Polyclonal Antibody detects endogenous levels of CLIP-170 protein.
Storage and Stability:	-20°C/1 year