

Antibodies Product Datasheet

Product Name:	CLIP-115 Antibody
Catalog Number:	C14916
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLIP2. AA range:997-1046
Concentration:	1 mg/ml
Gene Name:	CLIP2
Protein Name:	CAP-Gly domain-containing linker protein 2
Other Name:	CLIP2; CYLN2; KIAA0291; WBSCR3; WBSCR4; WSCR4; CAP-Gly domain-containing linker protein 2; Cytoplasmic linker protein 115; CLIP-115; Cytoplasmic linker protein 2; Williams-Beuren syndrome chromosomal region 3 protein; Williams-Beuren syndro
Human Swiss Prot No.:	Q9UDT6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLIP-115 Polyclonal Antibody detects endogenous levels of CLIP-115 protein.
Storage and Stability:	-20°C/1 year